


Made-in-The-USA

With Globally Sourced Materials

We Make Masks Right

3-Ply face mask (blue)



- High Quality Filtration
- Adjustable Nose Strip
- Welded Ear loops
- Breathable & Comfortable



Sponsor:
Scott Yin
Mid-America Safety Corp
800 Green Leaf Ave.
Elk Grove Village, IL 60007

Bacterial Filtration Efficiency (BFE) Final Report


Test Article: MB-a-1
Purchase Order: QT 0722-1
Study Number: 1326233-S01
Study Received Date: 30 Jul 2020
Testing Facility: Nelson Laboratories, LLC
6280 S. Redwood Rd.
Salt Lake City, UT 84123 U.S.A.
Test Procedure(s): Standard Test Protocol (STP) Number: STP0004 Rev 18
Deviation(s): None

Summary: The BFE test is performed to determine the filtration efficiency of test articles by comparing the bacterial control counts upstream of the test article to the bacterial counts downstream. A suspension

Results:

Test Article Number	Percent BFE (%)
1	>99.9
2	>99.9
3	>99.9 ^a
4	>99.9
5	>99.9 ^a

^a There were no detected colonies on any of the Andersen sampler plates for this test article.



50 pcs/per box; 20 boxes/carton;

